DC Bead^{MI}

embolic Drug-Eluting Bead
INSTRUCTION FOR USE

STERILE
SINGLE USE ONLY
NON-PYROGENIC

Sterilised by steam

Do not use if the package is opened or damaged

DESCRIPTION:

DC Bead^{MI} are precisely calibrated, hydrogel embolic drug-eluting beads. DC Bead^{MI}, produced from polyvinyl alcohol, are biocompatible, hydrophilic, and non resorbable.

DC Bead^{MI} are capable of loading and eluting irinotecan and are available in the following size:

Nominal Bead Size	Label Colour	Upon loading with irinotecan, DC Bead ^{MI} undergo a decrease in size, up to 30% when loading at 50mg/ml.
70 - 150µm	Black and Yellow	

PRESENTATION:

- 10ml glass vial.
- Each vial contains approximately 2ml of DC Bead^{MI} in non-pyrogenic, sterile, physiological buffered saline. Total volume of saline and DC Bead^{MI} is approximately 8ml.
- The vial is stopper sealed by an aluminium cap equipped with a colour-coded lid.
- Each vial is intended for single patient use only. Do not resterilise. Discard any unused material.

INDICATIONS:

DC Bead^{MI} is primarily intended as an embolic agent to treat vessels supplying malignant colorectal cancer metastasised to the liver (mCRC). DC Bead^{MI} is compatible with irinotecan, which can be loaded prior to embolisation and then, as a secondary action, elute a local, controlled and sustained dose to the mCRC after embolisation.

CONTRAINDICATIONS – DC Bead^{MI}:

- Patients intolerant to vascular occlusion procedures.
- Vascular anatomy that precludes catheter placement or emboli injection.
- Presence or likely onset of vasospasm.
- Presence or likely onset of haemorrhage.
- Presence of severe atheromatous disease.
- Presence of feeding arteries smaller than distal branches from which they emerge.
- Presence of patent extra-to-intracranial anastomoses or shunts.
- Presence of collateral vessel pathways potentially endangering normal territories during embolisation.
- Presence of end arteries leading directly to cranial nerves.
- Presence of arteries supplying the lesion not large enough to accept DC Bead^{MI}.
- Vascular resistance peripheral to the feeding arteries precluding passage of DC Bead^{MI} into the lesion.

Do not use DC Bead^{MI} in the following applications:

- Embolic agent in non-malignant tumours.
- Embolic agent in the presence of large diameter AV shunts.
- Embolic agent of arteriovenous shunts (ie. where the blood does not pass through the arterial/capillary/venous transition but directly from artery to vein).
- Any vasculature where DC Bead^{MI} embolic agent could pass directly into the internal carotid artery or other non-target territories.

CONTRAINDICATIONS – IRINOTECAN:

- See irinotecan package insert for contraindications regarding use.

WARNING: Studies have shown that DC Bead^{MI} do not form aggregates and, as a result, penetrate deeper into the vasculature as compared to similarly sized PVA particles.

CAUTIONS:

- Do not use if the vial or packaging appear damaged.
- DC Bead^{MI} is not intended for use with doxorubicin
- On addition of non-ionic contrast/water mixture to loaded beads some irinotecan will be eluted over time. If the beads are not used immediately, up to 10mg irinotecan may be present in the contrast/water mixture. If this occurs a small dose of irinotecan may be available systemically at time of delivery.
- Do not use irinotecan loaded beads with contrast agents containing salts (eg. Calcium chloride)
- Ensure that DC Bead^{MI} is an appropriate size for the intended vasculature.
- Consider upsizing to a larger size of DC Bead in the presence of AV shunts or if angiographic evidence of embolisation does not appear quickly during delivery.
- Consideration should be given to Tc99m-MAA scanning if there is suspicion of AV shunting.
- Monitor patients carefully for signs of non-target embolisation such as hypoxia or CNS changes.
- Embolization with DC Bead^{MI} should only be performed by a physician with appropriate interventional occlusion training in the region intended to be embolised.
- The maximum amount of irinotecan that can be loaded is 100mg irinotecan per 2ml vial of DC Bead^{MI}. Exceeding this amount may lead to some irinotecan remaining free in solution. This free solution should be removed prior to use to prevent the patient receiving the excess dose as a bolus.

POTENTIAL COMPLICATIONS:

- Undesirable reflux or passage of DC Bead^{MI} into normal arteries adjacent to the targeted lesion or through the lesion into other arteries or arterial beds.
- Non-target embolisation.
- Pulmonary embolisation.
- Ischaemia at an undesirable location.
- Capillary bed saturation and tissue damage.
- Ischaemic stroke or ischaemic infarction.
- Vessel or lesion rupture and haemorrhage.
- Neurological deficits including cranial nerve palsies.
- Vasospasm.
- Death.
- Recanalisation.
- Foreign body reactions necessitating medical intervention.
- Infection necessitating medical intervention.
- Clot formation at the tip of the catheter and subsequent dislodgement causing arterial thromboembolic sequelae.

DRUG LOADING INSTRUCTIONS:

DC Bead^{MI} is suitable for loading with irinotecan solution only.

A loading solution of 20mg/ml irinotecan is recommended.

To obtain a final loading of 100mg irinotecan per 2ml vial of DC Bead^{MI}:

- Remove as much saline as possible from a vial of DC Bead^{MI} using a syringe with a small gauge needle.
- Using a syringe and needle add the 5ml of irinotecan solution (20mg/ml) directly to the vial of DC Bead^{MI}.
- Agitate the DC Bead^{MI}/irinotecan solution gently to encourage mixing, then allow to stand. Agitate the solution at intervals. The beads will turn a turquoise colour as the loading progresses
- Loading will take a minimum of 2 hours
- Prior to use, transfer the DC Bead^{MI} loaded with irinotecan to a syringe, remove and discard any excess supernatant. Add 5ml water for injection (NOT sodium chloride) and 5 - 10ml of non-ionic contrast medium per ml of DC Bead^{MI}. Mix gently to give a good suspension. Deliver suspension at a rate of 1ml per minute.
- A dose of up to 100mg irinotecan per 2ml vial of DC Bead^{MI} can be loaded.
- Up to 4ml of beads loaded with irinotecan have been documented in clinical use (Fiorentini, In vivo 21, 1085, 2007).

STORAGE OF DC Bead^{MI}:

In order to minimise the risk of microbiological contamination DC Bead^{MI} should be prepared under controlled aseptic conditions. As the preparation and loading conditions of DC Bead^{MI} are outside of the manufacturers control, once the DC Bead^{MI} vial has been pierced, the allocation of a shelf life longer than 4 hours if used at room temperature or 24 hours if stored in a refrigerator at 2 - 8°C is the responsibility of the user.

DC Bead^{MI} loaded with irinotecan is physically and chemically stable for 14 days if stored in a refrigerator at 2 - 8°C. Once mixed with contrast media DC Bead^{MI} loaded with irinotecan must be used immediately.

DELIVERY INSTRUCTIONS:

- Carefully evaluate the vascular network associated with the lesion using high resolution imaging prior to beginning the embolisation procedure.
 - Ensure DC Bead^{M1} is suitable for the pathology to be treated (i.e. vascular target/vessel size).
 - Choose a delivery catheter based on the size of the target vessel. DC Bead^{M1} can tolerate temporary compression of 20% to 30% in order to facilitate passage through the delivery catheter.
 - Introduce the delivery catheter into the target vessel according to standard techniques.
 - DC Bead^{M1} is not radio-opaque. It is recommended to monitor the embolisation under fluoroscopic visualisation by adding the desired amount of contrast medium to the suspension fluid.
- i. Take care to ensure proper suspension of the DC Bead^{M1} in the contrast medium to enhance distribution during injection.
 - ii. Draw the DC Bead^{M1} into a syringe needle of a size greater than or equal to 19G (1.07mm).
 - iii. Slowly inject DC Bead^{M1} into the delivery catheter under fluoroscopic visualisation while observing the contrast flow rate. Exercise conservative judgment in determining the embolisation endpoint.
- Upon completion of the treatment, remove the catheter while maintaining gentle suction so as not to dislodge DC Bead^{M1} still within the catheter lumen.
 - Discard any unused DC Bead^{M1} loaded with irinotecan.

The following table shows the sizes of catheter which are compatible with delivering the product.

Microcatheter Size	Loading with Irinotecan	70 - 150µm
5Fr	Unloaded	Y
	Loaded	Y
3.0 Fr	Unloaded	Y
	Loaded	Y
2.7 Fr	Unloaded	Y
	Loaded	Y
2.4Fr	Unloaded	Y
	Loaded	Y
2.3 Fr	Unloaded	Y
	Loaded	Y*
1.8 Fr	Unloaded	Y
	Loaded	Y*

*Not tested but expected to deliver without any issues

CONSERVATION AND STORAGE:

Store unopened DC Bead^{M1} in a cool, dry and dark place in its original packaging.

Use by the date indicated on the vial label.

Do not freeze.

PACKAGE LABEL:

REF Catalogue number	 Steam Sterilised	 Protect from moisture
LOT Batch number/Lot number	 Use before/Expiry	 Attention see instructions for use
 Do not reuse	 Protect from light	 Do not freeze

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